## SECURITY HIGH NOBLE WHITE PORCELAIN ALLOY

Jensen Dental's **SECURITY** Alloy offers cost-conscious customers the ability to purchase a high-noble PFM dental alloy at an economical price. Since SECURITY has the smallest amount of gold possible to be classified as a "high noble" alloy (40%), it may be less expensive than other products that contain more gold. SECURITY's properties make it ideal for just about any PFM application. A strong alloy, it's great for single units, short- and long-span bridges.

PROPERTIES			
Melting Range	2120°-23	300°F (1160°-1260°C)	
Coefficient of Thermal Expansion			
from 25°C	to 500°C	14.2x10 <sup>-6</sup> K <sup>-1</sup>	
from 25°C	to 600°C	14.5x10 <sup>-6</sup> K <sup>-1</sup>	
Density		13.0 g/cm <sup>3</sup>	

40%		
40%		
9%		
6%		
3%		
2%		
Contains less than 1% Rhenium		
е		

AFTER PORCELAIN FIRING		
Hardness	255 HV	
Tensile Elongation	20%	
TensileYield Strength/Proof Stress	78,230 psi (540 Mpa)	
Ultimate Tensile Strength	110,230 psi (772 MPa)	
Modulus of Elasticity	17.4 x 10 <sup>6</sup> psi (120,000 MPa)	

## **PROCESSING TECHNIQUE**

WAXING 0.3mm for single units and 0.5mm for bridgework

**INVESTMENT** Phosphate-bonded high heat investment with or without carbon

BURNOUT 1550°F (850°C)

**MELTING AND** Wind the casting arm one more time than normal (3-4 winds). CASTING Multi-orifice torch: 10 psi fuel and 20 psi oxygen (0.7 bars fuel and 1.4 bars oxygen). Heat until cloudy surface clears before releasing the casting arm. DO NOT OVERHEAT. DO NOT USE CASTING FLUX. Casting temperature ~2400°F (1315°C).

FINISHING Finish with aluminum oxide stones. Blast porcelain receiving surfaces with non-recycled aluminum oxide. Clean in ultrasonic for 10 minutes.

## OXIDATION Oxidize from 1200°F (650°C) to 1850°F (1010°C) at 145°F/min (80°C/min) in vacuum with 5 minutes hold time. Bench cool. Proceed with normal opaque technique.

SOLDERS AND FLUX Pre-solder: PWS or LX Solder Post-solder: 1400 Solder Flux: Brown Fluoride Flux (for both pre & post soldering)

> 5225Y r2



Jensen Industries, Inc. 50 Stillman Road North Haven, CT 06473 USA EC REP

Jensen GmbH Gustav – Werner Str. 1 D-72555 Metzingen Freecall: 0800-8 57 32 30 info@jensendental.de www.jensendental.de

Jensen Technical Service is available 8:00 AM - 8:00 PM EST by calling 1-800-243-2000

