



PMMA Resin Material

Mono | Multi | Gum | Transparent |



Multi-layer Temporary Restoration PMMA

CAD/CAM Anterior esthetic temporary restoration material, with revolutionary shade gradience technology, demarcation of shades realization, Multi-layer PMMA satisfies shade gradience from incisal to bottom. Excellent performance, stability, translucency and good bio-compatibility, multi-layer PMMA delivers no-stimulus features to the oral tissues.



Personalized high esthetics

Protect abutment from stimulus, polishing friendly, high gloss

Shade stability, and anti-shade fade



Indication:



- Temporary Crown | • Temporary Bridge | • Temporary Full Arch | • Full Contour Screw Retained bridge

Specification:

Bending Strength	≥120MPa	Dimension	12/14/16/18/20/22/25/30
System	Wieland/Amann/Zirkon/Sirona	Shade	Vita 16 shades/OM1/OM2/OM3
Water Sorption	≤25[μg/mm ³]	Solubility	≤5[μg/mm ³]

Aidite / PMMA Resin Material