


KATANA™ Zirconia



STML Super Translucent Multi-Layered. Ideal for up to 3 units posterior bridges with a well balanced combination of chromatic and gradational translucency, which reproduces aesthetic enamel and dentin effects.

HTML High Translucent Multi-Layered. Suitable for anterior and/or posterior monolithic restorations that need high strengths and where the effect of the abutment color needs to be suppressed to a minimum. It can also be used for single crown restorations and for long-span bridges. Furthermore, it is an excellent framework material when used in combination with CerabienZR porcelain.

Mechanical Property



Date source: Kuraray Noritake Dental Inc. The numerical value varies depending on the conditions.

KATANA™ ZIRCONIA

KATANA™ ZIRCONIA is the trade mark of high performance zirconia discs for crowns, bridges and frames with highly innovative features.

With the KATANA™ ZIRCONIA discs you have excellent machinability, which avoids dipping during your production process. These zirconias contain advanced ceramic technology, which results in longevity. They show aesthetic translucency, exceptional shrinkage consistency on every distribution point within the disc, as well as high flexural strength.

