

Temp Esthetic

Multi-layered composite PMMA for long term temporaries in the esthetic zone.

Finally. A provisional solution for the anterior esthetic zone.

A highly esthetic, long term provisional solution offering unparalleled shade gradation which accurately represents the characteristics of natural teeth.

Homogeneous shade gradation.

Composed of 17 layers of incisal gradient, Temp Esthetic consists of a hybrid composite PMMA structure made of acrylic resins and organically modified ceramics. The result is a harmonious distribution of light, opalescence and fluorescence.

